

Leneta Company, Inc.

15 Whitney Rd.

Mahwah, NJ, 07430 USA

Phone: 201.847.9300 • **Fax:** 201.848.8833**Email:** sales@leneta.com • **Website:** www.leneta.com**Form 2DX Brushout Card**

Designed for informal brushout applications, the paper stock is thicker than regular chart paper to give greater rigidity for more convenient handling. They are also used widely for drawdowns and colorimetric measurements.

Form 2DX is a black and white brushout card. The top is black and the bottom is white with overall dimensions of 3-7/8 x 6 in (98 x 152 mm).

Specifications

Color	Black (Top) & White (Bottom)
Size	3-7/8 x 6 in 98 x 152 mm
Box Quantity	500
Boxes per Case	4
Box Weight	7 lb
Nominal Thickness	20 mil 0.5 mm