

Leneta Company, Inc.

15 Whitney Rd.

Mahwah, NJ, 07430 USA

Phone: 201.847.9300 • **Fax:** 201.848.8833**Email:** sales@leneta.com • **Website:** www.leneta.com**Form 3B Opacity Chart**

Form 3B is a black and white sealed opacity chart. The top and bottom are white, middle section is black. Overall dimensions are 7-5/8 x 11-3/8 in (194 x 289 mm).

The header includes a writable area for recording:

Date, Time
Formula
Batch Number
Notes

A 1/4" diameter hole is centered just under the top of the chart.

Specifications

Color	Black & White
Size	7-5/8 x 11-3/8 in 194 x 289 mm
Box Quantity	250
Boxes per Case	4
Box Weight	9 lb
Minimum Order	1 Box