


Leneta Company, Inc.
 15 Whitney Rd.
 Mahwah, NJ, 07430 USA
Phone: 201.847.9300 • **Fax:** 201.848.8833
Email: sales@leneta.com • **Website:** www.leneta.com



Form 3NT-1, Printing Ink Drawdown Sheets

Vellum Opaque
 This is a neutral white, vellum opaque paper, unsealed. **Contains no optical brighteners.**

Specifications	
Sheet Size	5 x 7-5/8 in 127 x 194 mm
Paper	Unwatermarked, non-fluorescent
Ink	Jet black, Non-bleeding
Paper Type	Vellum Opaque
Shade	Neutral White
Basis Weight ¹	60 lb lb
Grammage	89 g/m²
Density	0.70 g/cm³
Boxes per Case	4
Box Weight	6 lb lb
Minimum Order	1 box

¹ Ream of 500 basis sheets.