

Leneta Company, Inc.
15 Whitney Rd.
Mahwah, NJ, 07430 USA
Phone: 201.847.9300 • **Fax:** 201.848.8833
Email: sales@leneta.com • **Website:** www.leneta.com



Form 9B Opacity-Display Chart

Form 9B is a black and white sealed opacity-display chart. The top most area is strip, then black and then white with overall dimensions of 7-5/8 x 11-3/8 in (194 x 289 mm).

Specifications

Color	Black & White
Size	7-5/8 x 11-3/8 in 194 x 289 mm
Box Quantity	250
Boxes per Case	4
Box Weight	9 lb