Dial Reading Viscometer

...our original now with electronic drive

The Worldwide Standard Viscometer

Easy-to-Select Speeds

Analog Display
-Shows % Torque
-Use Factor Finder to convert reading
to centipoise

Simple-to-Use for easy setup

Available in Explosion Proof
U.L. Class 1, Group D lodcations
(without Electronic Drive)

Accuracy: ±1.0%
-Repeatability: ±0.2%

2-Year Limited Warranty

Accuracy: ±1.0% of range –Repeatability: ±0.2%



BROOKFIELD **V**ISCOMETERS

United States: T: 800.628.8139 or 508.946.6200 F: 508.946.6262 www.brookfieldengineering.com United Kingdom: T: 01279 451774 F: 01279 451775 www.brookfield.co.uk

What's Included?

Instrument

6 spindles (RV/HA/HB)

or 4 spindles (LV)

Spindle Guard Leg*

Lab Stand (Model A)

Carrying Case

Optional Accessories

Viscosity Standards

RV/HA/HB-1 Spindle

Quick Action Lab Stand

Temperature Bath

Small Sample Adapter

UL Adapter

Thermosel

Helipath Stand with T-bar Spindles

Spiral Adapter

DIN Adapter

Quick Connect/Extension Links

Vane Spindles



Easy Speed Adjustment and On/Off Control

Dial Viscometer

Model Selection Table Brookfield Dial Viscometers

VISCOSITY RANGE cP(mPa•s)

SPEEDS

MODEL	Min.	Max.	RPM	Number of Increments
LVT	1 [†]	2M	.3-60	8
RVT	100 ^{††}	8M	.5-100	10
HAT	200 ^{††}	16M	.5-100	10
HBT	800 ^{††}	64M	.5-100	10

† 1 cP achieved with UL Adapter accessory. 15 cP on LV with standard spindles.

†† Minimum viscosity is achieved with optional RV/HA/HB-1 spindle.

M = 1 million cP = Centipoise mPa•s = Millipascal/seconds

Applications

Consider application and viscosity range when selecting model (LV, RV, HA, HB)

LV SERIES – LOW VISCOSITY

Adhesives (solvent base) Biological Fluids Chemicals Dairy Products Inks

Juices
Oils / Solvents
Pharmaceuticals
Polymer Solutions

RV SERIES – MEDIUM VISCOSITY

Adhesives (hot melt) Asphalt (SHRP) Ceramic Slurries Food Products Gums

Ointments / Creams Paints / Inks / Coatings Paper Pulp

Shampoos / Lotions

HA/HB SERIES – HIGH VISCOSITY

Caulking Compounds / Putty

Chocolate Epoxies

Inks (ballpoint, offset, lithographic)

Pastes
Peanut Butter
Roofing Compounds

Sealants Tars



11 Commerce Boulevard Middleboro, MA 02346-1031 USA T: 800.628.8139 or 508.946.6200 F: 508.946.6262 W: www.brookfieldengineering.com

^{*}Not applicable to HA or HB versions