

KU-3™ Viscometer

for Paints, Coatings and Inks

ASTM D562 Compatible
(industry specification)

NEW Magnetic Spindle Coupling

NEW User Interface

**NEW Single Piece Pint and 1/2 Pint
Can Adapter**

Easy to use

no weights, simplifies an
established test procedure

LED Display Info:

- Krebs Units
- Gram Units
(Weight)
- Centipoise*

Select Krebs or Grams or Centipoise

Lock-In Test Results
with hold button

Accuracy: $\pm 1.0\%$ of range

Repeatability: $\pm 0.5\%$

Standard Krebs Spindle

Measurement range:

40 to 141 KU, 32 to 1099 gm, and 27
to 5274 cP*

Adapter provided

for 1/2 Pint, Pint,
& Quart Containers



What's Included?

Instrument
Krebs-type Spindle (KU-1030)
Adapter for Quart, Pint and
Half-Pint Can

Optional Accessories

Paste Spindle (KU-1031)

Applications

Paints
Coatings
Adhesives
Inks
Pastes

*Centipoise values based on the conversion from Krebs Units
as defined in ASTM D562.

ISO 9001:2008 CERTIFIED

BROOKFIELD
AMETEK®