

Vision on quality www.tqc.eu



TQC Scrub Abrasion and Washability Tester

AB6000, AB6010

Used to test the resistence of paint, vanish or coatings to scratching, wearing, and color loss due to wet or dry abrasion, by simulating everyday wear from cleaning actions or general use. The TQC Scrub Abrasion and Washability Tester is suitable for a wide range of scrub, abrasion and washability tests that are in use.





Flexible adaptation to many tests

The TQC Scrub Abrasion and Washability Tester allows for the most flexible adaptation of the wide range of scrub, abrasion and washability tests that are in use. First the base apparatus has to be selected from two models. The TQC Scrub Abrasion and Washability Tester - Basic (Art. Nr. AB6010) without pumps for dry tests, and the TQC Scrub Abrasion and Washability Tester - Standard (Art. Nr. AB6000) with pumps for dry and wet tests. Each base unit is completely flexible and allows for easy exchange of the tools used for testing.

A variety of tests are possible, with just small modifications to make the device suitable for other tests. This modular system significantly decreases the costs for laboratories to comply to many standards.

Applications



Detergents



Abrasives



Seat Covers



Phone Screens

and many more...

Features



Preset range option



Full color display



Triple i Control (Intelligent Illumination Interface)



Easy to update software



Manual movement



Compatible with a large range of accessories

