

Leneta Company, Inc.

15 Whitney Rd.

Mahwah, NJ, 07430 USA

Phone: 201.847.9300 • **Fax:** 201.848.8833**Email:** sales@leneta.com • **Website:** www.leneta.com**Plain White Charts**

Special sizes available upon request.

Form WH is a plain white chart with overall dimensions of 11-1/4 x 17-1/4 in (286 x 438 mm).

Specifications

Color	Plain White
Size	11-1/4 x 17-1/4 in 286 x 438 mm
Box Quantity	125
Boxes per Case	4
Weight per Box	11 lb
Nominal Thickness	14 mil 0.36 mm